



Did you know that ChemInstruments offers calibration and upgrade on all of our hand rollers?

Hand rollers help to lay the foundation of accurate product testing
[pictured: HR-100]



Why do I need to have my product calibrated? You may be wondering why you need to have something you've purchased from ChemInstruments calibrated. Over time, wear and tear from testing take their toll and machines and accessories will no longer be as accurate as when first purchased. We test for changes in weight, width, and hardness of the roller.

Why choose calibration? Calibration is a cost effective alternative to replacing your products annually. Additionally, the certificate you receive is good for another year of cooperating with ASTM and PSTC standards.

How often should I do this? We recommend calibrating your items once a year for optimal performance.



[pictured: HR-FINAT]

Certificates and customization available upon request



PSTC
PRESSURE SENSITIVE TAPE COUNCIL
A bond you can trust.



SERVING
INDUSTRY
SINCE 1986

510 Commercial Drive
Fairfield, OH 45014
p: 513.860.1598
f: 513.860.1597
sales@cheminstruments.com
www.cheminstruments.com